

DASSIET

UCAST[®]

UCAST is a pre-made super-splint that consists of two parts, a 100% biodegradable splint and a cohesive fabric bandage. Apply in 5 minutes, remold up to 25 times and recycle fully.



U-WRIST



U-THUMB-SH



U-HUM-L/R



U-MC



U-THUMB-L



U-ANKLE

PRODUCT		SIZE	BATCH
U-WRIST	Wrist splint (dorsal/volar)	One size fits most	10 pcs
U-THUMB-SH	Short thumb splint	One size fits most	10 pcs
U-HUM-L	Humerus orthosis, left	One size fits most	10 pcs
U-HUM-R	Humerus orthosis, right	One size fits most	10 pcs
U-MC	Metacarpal splint	One size fits most	10 pcs
U-THUMB-L	Long thumb splint	One size fits most	10 pcs
U-ANKLE	Ankle brace	One size fits most	10 pcs

UCAST[®]



U-FINGER



U-LONGF



U-MALLET



U-PATELLA



U-KNEE

PRODUCT		SIZE	BATCH
U-FINGER-SM	Finger buddy loop	Small-Medium	10 pcs
U-FINGER-ML	Finger buddy loop	Medium-Large	10 pcs
U-LONGF	Long finger splint	One size fits most	4 pcs
U-MALLET	Mallet finger splint (dorsal/volar)	One size fits most	10 pcs
U-PATELLA	Patella support	One size fits most	1 pc
U-KNEE	Long knee splint	One size fits most	1 pc
U-TOE	Big toe / hallux valgus splint	One size fits most	10 pcs

WOODCAST™

Woodcast is a heat-moldable biomaterial ideal for making casts, splints and orthoses. It's safe, X-ray transparent, remoldable and selfattaching. Use our heating devices to work with Woodcast.

WOODCAST 1 MM RIBBON	REF	FINISH	SIZE	BATCH	COLOR
For finger and thumb splinting and reinforcement.	D4500	Solid	1mm x 4cm x 5m	2 pcs	
	D8750	Solid	1mm x 8cm x 5m	1 pcs	



WOODCAST 2MM	REF	FINISH	SIZE	BATCH	COLOR
For all upper limb applications and smaller lower limb applications.	C1434	Vented	2mm x 14.5cm x 34cm	20 pcs	
	C1480	Vented	2mm x 14.5cm x 80cm	10 pcs	



WOODCAST 2MM SOFT	REF	FINISH	SIZE	BATCH	COLOR
Ideal for hand therapy and focused rigidity casting.	F1434	Vented	2mm x 14.5cm x 34cm	20 pcs	
	F1480	Vented	2mm x 14.5cm x 80cm	10 pcs	



WOODCAST 2 MM FOLDED	REF	FINISH	SIZE	BATCH	COLOR
Folded version of 80cm long Woodcast 2mm that fits in the Express Heater.	C1480F	Vented	2mm x 14.5cm x 80cm	2 pcs	



WOODCAST™

WOODCAST 2MM 4IN1	REF	FINISH	SIZE	BATCH	COLOR
Pre-cut splint for most common upper limb applications.	C1441	Vented	2mm x 14,5cm x 29cm (pre-cut in shape)	20 pcs	
					

WOODCAST ONESTEP	REF	FINISH	SIZE	BATCH	COLOR
Pre-cut splint for ankle back-slabs and above elbow casts.	C1465 (Elbow)	Vented	2mm x 14.5cm x 65cm	10 pcs	
	C1565 (Ankle)	Vented	4mm x 15cm x 65cm	5 pcs	
					

WOODCAST 4MM	REF	FINISH	SIZE	BATCH	COLOR
Thick splint that provides good support. Used for lower limb applications.	B1040	Vented	4mm x 10cm x 40cm	20 pcs	
	B1080	Vented	4mm x 10cm x 40cm	10 pcs	
	B1240	Vented	4mm x 12.5cm x 40cm	20 pcs	
	B1280	Vented	4mm x 12.5cm x 80cm	10 pcs	
	B1540	Vented	4mm x 15m x 40cm	20 pcs	
	B1580	Vented	4mm x 15m x 80cm	10 pcs	
					

WOODCAST 4MM FOLDED	REF	FINISH	SIZE	BATCH	COLOR
Pre-cut splint for most common upper limb applications.	B1580F	Vented	4mm x 15m x 80cm	8 pcs	
					

WOODCAST™

ANKLE BRACE	REF	SIZE	BATCH	COLOR
Moldable ankle brace, one size fits all.	YL2808	2mm x 14.5cm x 80cm	2 pcs	Black




WOODCAST TCC-KIT	REF	SIZE	BATCH	COLOR
One kit that contains everything needed for one Total Contact cast. Available in two thicknesses.	TCC6020	6mm thickness	1 kit	
	TCC8020	8mm thickness	1 kit	



PROCAST™

PROCAST is our line of field-tested accessories for safe and effective treatment, designed to work perfectly with WOODCAST and UCAST.

STOCKINETTE	REF	SIZE	BATCH	COLOR
<p>Soft and moisture-wicking stockinette for protecting the skin.</p> 	X2522	2.5cm x 22.8m	1 pc	White
	X5022	5cm x 22.8m	1 pc	White
	X7522	7.5cm x 22.8m	1 pc	White
	X1022	10cm x 22.8m	1 pc	White
	X1522	15cm x 22.8m	1 pc	White

ADHESIVE AIRPADDING	REF	SIZE	BATCH	COLOR
<p>Our unique padding material, that is extremely lightweight. Breathable and easy to work with.</p> 	X4350A	40cm x 3.5m	1 pc	Black

EDGING TAPE	REF	SIZE	BATCH	COLOR
<p>For finishing any hard edges.</p> 	X0315	3cm x 7.5m	2 pcs	Black

HOOK	REF	SIZE	BATCH	COLOR
<p>Adhesive Hook for orthoses that offers a unique, strong hold.</p> 	X5250	50mm x 22.8m	1 pc	Black

PROCAST™

ELASTIC LOOP	REF	SIZE	BATCH	COLOR
Elastic loop for making orthoses. Use together with Hook.	X0250	25mm x 7.5m	2 pcs	Black
	X5252	50mm x 16m	1 pc	Black



LOOP	REF	SIZE	BATCH	COLOR
Loop for making orthoses. Use together with Hook.	X5251	50mm x 22.8m	1 pc	Black



SOFT SLING	REF	SIZE	BATCH	COLOR
Supports the arm without chafing the neck. Comes with 28 Y-hooks.	X5524	2 x 12m	2 pcs	Black



SCISSORS	REF	PRODUCT	BATCH	REMARKS
Durable scissors that stay sharp for a long time.	MPX2020	Economy scissor	1 pc	
	MPX1980	TC Lister Scissor	1 pc	Can be sharpened
	MP1811	Curved Nibbler Scissor	1 pc	
	X1811A	Curved Nibbler Scissor TC	1 pc	Can be sharpened



PROCAST™

EXPRESS HEATER

Faster heating device for Woodcast and UCAST.



REF

A002

SIZE

65cm x 23cm x 39cm

BATCH

1 pc

COLOR

White

FLATBED HEATER

Large flatbed heater for Woodcast and UCAST.

REF

A003

SIZE

10cm x 44cm x 88cm

BATCH

1pc

COLOR

Grey

FREEZE SPRAY

Used to cool down splint faster.



REF

X0400

SIZE

400 ml

BATCH

1pc

COLOR

-